

Sample Name :[D-SER]33-Mazdutide  
Cat No. :RP30899  
Lot No. :E9748225H  
Time Processed :12:57:16  
Month-Day-Year Processed :08/30/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
Total Flow:1 ml/min  
Wavelength:220 nm  
<<LC Time Program>>


| Time  | Module     | Command | Value |
|-------|------------|---------|-------|
| 0.01  | Pumps      | B.Conc  | 5     |
| 25.00 | Pumps      | B.Conc  | 65    |
| 25.01 | Pumps      | B.Conc  | 95    |
| 27.00 | Pumps      | B.Conc  | 95    |
| 27.01 | Pumps      | B.Conc  | 5     |
| 35.00 | Pumps      | B.Conc  | 5     |
| 35.01 | Controller | Stop    |       |

<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm  
Equipment: ZJ19010141

<Chromatogram>



<Peak Table>

| Detector A 220nm |           |         |        |         |
|------------------|-----------|---------|--------|---------|
| Peak#            | Ret. Time | Area    | Height | Area%   |
| 1                | 21.146    | 65649   | 7572   | 1.309   |
| 2                | 21.332    | 62945   | 6748   | 1.255   |
| 3                | 21.743    | 4878766 | 463114 | 97.256  |
| 4                | 23.114    | 9082    | 1057   | 0.181   |
| Total            |           | 5016442 | 478492 | 100.000 |