

Sample Name :[D-GLU]16-Mazdutide

Cat No. :RP30900

Lot No. :E9749525H

Time Processed :17:11:52

Month-Day-Year Processed :09/24/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010325

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.467	37294	871	0.217
2	14.412	23484	815	0.137
3	16.272	10602	1038	0.062
4	16.550	33039	2156	0.192
5	16.852	37041	3177	0.216
6	17.188	130066	8066	0.757
7	17.735	16862958	1291069	98.132
8	18.983	10330	1640	0.060
9	19.315	39112	1767	0.228
Total		17183927	1310598	100.000