

Sample Name :[D-GLU]16-Mazdutide
Cat No. :RP30900
Lot No. :E9750325H
Time Processed :12:21:50
Month-Day-Year Processed :08/30/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm
<<LC Time Program>>


Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010141

<Chromatogram>



<Peak Table>

Detector A 220nm				
Peak#	Ret. Time	Area	Height	Area%
1	8.318	8835	238	0.127
2	13.176	9401	687	0.135
3	20.407	31559	2095	0.453
4	20.851	66103	4674	0.949
5	21.451	6749486	639713	96.940
6	22.065	97151	5531	1.395
Total		6962536	652938	100.000