


Sample Name :[D-ALA]19-Mazdutide
Cat No. :RP30908
Lot No. :E9765225H
Time Processed :18:36:24
Month-Dav-Year Processed :09/23/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ17010504

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.617	40550	3683	0.704
2	19.850	20109	4201	0.349
3	20.217	5631818	561322	97.804
4	20.625	19125	1403	0.332
5	21.292	46647	4064	0.810
Total		5758249	574673	100.000