


Sample Name :[D-ALA]19-Mazdutide
Cat No. :RP30908
Lot No. :E9766125H
Time Processed :12:54:51
Month-Day-Year Processed :09/05/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: GK12010012

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm				
Peak#	Ret. Time	Area	Height	Area%
1	20.850	93320	8476	0.834
2	21.050	48595	5411	0.434
3	21.450	11052909	933149	98.732
Total		11194823	947036	100.000