


Sample Name :KRAS G12V Peptide (VVVGAVGVGK)
Cat No. :RP30229
Lot No. :E9766825K
Time Processed :23:57:56
Month-Day-Year Processed :11/18/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ24010222

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.325	24457	3800	1.285
2	9.117	1870580	257400	98.255
3	9.358	8769	1606	0.461
Total		1903805	262806	100.000