


Mass Spectrum



Sample Information		Interface	:ESI	Equipment	:ZJ22010150
Month-Day Processed :	11/18/25	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Time Processed :	21:14:57	CDL Temp	:250	Drying Gas Flow	:5 L/min
Injection Volume :	0.2	Block Temp	:200	T.Flow	:0.2 ml/min
Sample Name :	KRAS G12V Peptide (VVVGAVGVGK)				
Cat No. :RP30229	B.conc				
Lot No. :E9766825K					
Theoretical MW : 884.08					
Observed MW : 884.2					