


Mass Spectrum



Sample Information

Month-Day Processed : 08/28/25
 Time Processed : 8:43:54
 Injection Volume : 0.3
 Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. : RP11110
 Lot No. : E9790025H
 Theoretical MW : 2197.60
 Observed MW : 2197.2

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Equipment : ZJ21010035
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH