

Sample Name :TAT peptide
Cat No. :RP20256
Lot No. :E9795425H
Time Processed :9:56:16
Month-Day-Year Processed :08/28/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	0
5.00	Pumps	Solvent B Conc.	0
20.00	Pumps	Solvent B Conc.	30
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
30.00	Pumps	Solvent B Conc.	95
30.01	Pumps	Solvent B Conc.	0
38.00	Pumps	Solvent B Conc.	0
38.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.300	70361	42058	1.075
2	11.408	6248718	1067745	95.503
3	11.542	126368	56048	1.931
4	11.767	90357	13270	1.381
5	11.900	7163	1130	0.109
Total		6542967	1180251	100.000