


Mass Spectrum



Sample Information

Month-Day Processed : 08/27/25
 Time Processed : 22:20:06
 Injection Volume : 0.3
 Sample Name : TAT peptide
 Cat No. : RP20256
 Lot No. : E9795425H
 Theoretical MW : 1621.91
 Observed MW : 1621.2

Interface : ESI

Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : ZJ21010035

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H₂O/50%MeOH