


Mass Spectrum



Sample Information

Month-Day Processed :	03/31/23	Equipment	:ZJ22010150
Time Processed :	7:44:14 AM	Interface Bias	: +4.5 kV
Injection Volume :	0.1	Drying Gas Flow	:5 L/min
Sample Name :	RGD	T.Flow	:0.2 ml/min
Cat No. :	RP20297	B.conc	:50%H2O/50%MeOH
Lot No. :	E9856823C		
Theoretical MW : 1025.06			
Observed MW : 1025.2			