

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. :E9975026H
Time Processed :6:54:08
Month-Day-Year Processed :09/01/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.775	5665962	784753	95.286
2	6.967	255348	40170	4.294
3	7.133	7605	2029	0.128
4	7.433	17348	1972	0.292
Total		5946263	828924	100.000