

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E9976823B

Time Processed :4:39:00

Month-Day-Year Processed :02/14/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.828	9279	1770	0.089
2	9.086	195750	23001	1.867
3	9.543	37186	5345	0.355
4	9.761	10031027	1052551	95.697
5	10.013	147385	26431	1.406
6	10.302	53688	7388	0.512
7	10.482	7804	2774	0.074
Total		10482118	1119259	100.000