

## CERTIFICATE OF ANALYSIS

### Product Information

---

Product Name:	Broad Multi Color Pre-Stained Protein Standard
Cat. No.:	M00624-1250
Lot No.:	C70231907
Size:	250µl*5vial
Storage Condition:	-20°C

---

### Quality Control

---

Test Items	Specifications	Results
Physical Appearance	Clear liquid	Conformed
SDS-PAGE Testing	The bands are clear	Conformed
Protein Transfer Testing	The bands are clear	Conformed

---

### Caution:

For laboratory use only. Not for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 07/03/2019  
Product Department