

CERTIFICATE OF ANALYSIS

Product Information

Product Name:	Broad Multi Color Pre-Stained Protein Standard
Cat. No.:	M00624-250
Lot No.:	C70031904
Size:	250µl
Storage Condition:	-20 °C

Quality Control

Test Items	Specifications	Results
Physical Appearance	Clear liquid	Conformed
SDS-PAGE Testing	The bands are clear	Conformed
Protein Transfer Testing	The bands are clear	Conformed

Caution:

For laboratory use only. Not for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by:  Date: 04/03/2019
Product Department