

CERTIFICATE OF ANALYSIS

Product Information

Product Name:	Broad Multi Color Pre-Stained Protein Standard
Cat. No.:	M00624
Lot No.:	C70251910
Size:	250µl/1250µl
Storage Condition:	-20°C

Quality Control

Test Items	Specifications	Results
Physical Appearance	Clear liquid	Conformed
SDS-PAGE Testing	The bands are clear	Conformed
Protein Transfer Testing	The bands are clear	Conformed

Caution:

For laboratory use only. Not for household use. If you have any questions about the Certificate of Analysis, please contact our customer service representative at 1-877-436-7274 (Toll-Free), or 1-732-885-9188.

Certified by: *Xianli Yang* Date: 10/28/2019
Product Department