

**DATASHEET**

Version: 2016-04-19

**Cat. No.:** M01059**Size:** 200 µl**Category:** Immunology Reagents: cAMP-HRP**Storage:**

The conjugate can be stored for up to two years at -20°C or below.

**Applications:**

Working concentrations for specific applications should

be determined by the investigator. The appropriate concentrations may be affected by antibody affinity, the sensitivity of the method of detection, temperature, the length of the incubations, and other factors. The suitability of this product for applications other than those listed below has not been determined. The following concentration ranges are recommended starting points for this product.

ELISA: 1:10,000 - 1:20,000

Other: user-optimized