

**Pressinoic Acid****Cat. No.:** RP10408**Size:** 1 mg**Description:**

Pressinoic Acid is a peptide with potent corticotrophin-releasing activity. It is also an oxytocin inhibitor; it induces maternal behavior.

**Cas No:** 35748-51-7**Sequence (one-letter code):**

CYFQNC

**Sequence (three-letter code):**

{CYS}{TYR}{PHE}{GLN}{ASN}{CYS}

**Formula:** C<sub>33</sub>H<sub>42</sub>N<sub>8</sub>O<sub>10</sub>S<sub>2</sub>**Molecular Weight:** 774.86**Purity:** > 95%**Storage:**

Store the peptide at -20°C.

**Note:**

\*For Non-Clinical Research Use Only \*

