

Phe-Met-Arg-Phe**Cat. No.:** RP10814**Size:** 5 mg**Description:**

none

Sequence (one-letter code):

FMRF

Sequence (three-letter code):

{PHE}{MET}{ARG}{PHE}

Solubility: Soluble in water. The contents of this vial have been accurately determined. Both the stopper and the vial have been siliconized. Do not attempt to weight out a smaller portion of the contents.

Formula: C₂₉H₄₁N₇O₅S₁**Molecular Weight:** 599.8**Purity:** > 95%**Storage:**

Before using, store the peptide in the DRY form at 0-5°C. For best and repeatable results, rehydrate the peptide immediately before using. Do not re-freeze any unused portions.

Note:

*For Non-Clinical Research Use Only *

