

Rev03
Update: Jun,09,2025

DATASHEET

Exendin-4

Cat. No.: RP10874

Overview

Synonyms	Exendin4exenatide; Byetta
Description	Exendin-4 is a 39-amino acid peptide amide. Exendin-4, like Exendin-3, stimulates an increase in acinar cAMP without stimulating the release of amylase. Exendin-4 is a long-acting potent agonist of the glucagon-like peptide 1 (GLP-1). Exendin-4 is the active component of Byetta (exenatide) injection, which may improve glycemic control in people with type 2 diabetes mellitus and has the potential to reduce plasma glucose at least partly by a delay in gastric emptying as well as by reducing calorie intake. Exendin-4 enhances glucose-dependent insulin secretion by the pancreatic β -cell and suppresses inappropriately elevated glucagon secretion.
Cas No	141758-74-9
Sequence	{HIS}{GLY}{GLU}{GLY}{THR}{PHE}{THR}{SER}{ASP}{LEU}{SER}{LYS}{GLN}{MET}{GLU}{GLU}{GLU}{ALA}{VAL}{ARG}{LEU}{PHE}{ILE}{GLU}{TRP}{LEU}{LYS}{ASN}{GLY}{GLY}{PRO}{SER}{SER}{GLY}{ALA}{PRO}{PRO}{PRO}{SER}
Sequence Shortening	HGEGTFTSDLKQMEEEAVRLFIEWLKNGGPSSGAPPPS
Molecular Formula	C ₁₈₄ H ₂₈₂ N ₅₀ O ₆₀ S
C Terminal	Amidation
Molecular Weight	4186.66

Properties

Purity	> 95%
Solubility	Soluble in Ultrapure water under 1mg/ml
Form	Lyophilized
Storage	Store the peptide at -20°C.
Note	This peptide interacts with exendin receptor to increase pancreatic acinar cAMP. It has no secretagogue activity and does not bind to VIP receptors.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China