

## HA Peptide

**Cat. No.:** RP11735-10

**Size:** 10 mg

**Alias:** HA1 fragment 98-106Influenza Hemagglutinin (HA)  
Peptide

**Description:**

HA Peptide is useful for the displacement of HA-tagged fusion proteins from HA antibody in immunoassays. The successful inhibition of antibody binding by HA peptide demonstrates that the binding is specific. HA peptide recognized by anti HA; belongs to influenza hemagglutinin (HA) family and is responsible for attaching the virus to cell receptors and for initiating infection.

**Cas No:** 92000-76-5

**Sequence (one-letter code):**

YPYDVPDYA

**Sequence (three-letter code):**

{TYR}{PRO}{TYR}{ASP}{VAL}{PRO}{ASP}{TYR}{ALA}

**Formula:** C<sub>53</sub>H<sub>67</sub>N<sub>9</sub>O<sub>17</sub>

**Molecular Weight:** 1,102.15

**Purity:** > 95%

**Storage:**

Store at -20°C. Keep container tightly closed. Store in a cool dry place.

**Note:**

\*For Non-Clinical Research Use Only \*