

Rev01
 Update: Dec,28,2022

DATASHEET

CEF20

Cat. No.: RP11781

Overview

Synonyms	CEF20, Cytomegalovirus, CMV pp65 (495-503)
Description	CEF20, a CEF control peptide pool, is an HLA-A*0201-restricted epitope from cytomegalovirus pp65 (495-503). Such epitopes are useful in applications such as ELISPOT, intracellular cytokine, and CTL assays, and they also provide as net weight based on peptide content.
Sequence	{ASN}{LEU}{VAL}{PRO}{MET}{VAL}{ALA}{THR}{VAL}
Sequence Shortening	NLVPMVATV
Molecular Formula	C ₄₂ H ₇₄ N ₁₀ O ₁₂ S ₁
Molecular Weight	943.18

Properties

Purity	95%
Solubility	Soluble in water
Form	Lyophilized
Storage	Store at -20°C

Examples

