

PA (224-233), Influenza**DATASHEET**

Version: 2016-08-18

Cat. No.: RP19991**Size:** 0.5 mg**Description:**

This is a murine H-2 Db- and Kb-restricted immunodominant CTL epitope of influenza PR8 virus. The CD8+ T cells recovered by the bronchoalveolar lavage (BAL) from PR8-infected animals responded dominantly to stimulation with NP366-374 and PA224-233.

Sequence (one-letter code):

SSLENFRAYV

Sequence (three-letter code):

{Ser}{Ser}{Leu}{Glu}{Asn}{Phe}{Arg}{Ala}{Tyr}{Val}

Solubility: Soluble in water**Formula:** C₅₃H₈₀N₁₄O₁₇**Molecular Weight:** 1,185.29**Purity:** 95%**Storage:**

Store at -20°C

Note:

*For Non-Clinical Research Use Only *

