

Rev03
Update: Feb,10,2022

DATASHEET

NP (396-404)

Cat. No.: RP20090

Overview

Description	GenScript NP (396-404) is an immunodominant H-2Db-restricted LCMV peptide, restricted by H-2Db major histocompatibility complex class I molecules and presented to cytotoxic T.
Sequence	{Phe}{Gln}{Pro}{Gln}{Asn}{Gly}{Gln}{Phe}{Ile}
Sequence Shortening	FQPQNGQFI
Molecular Formula	C ₅₀ H ₇₁ N ₁₃ O ₁₄
Molecular Weight	1078.18

Properties

Purity	95%
Form	Lyophilized
Storage	Store GenScript NP (396-404) at -20°C.