

**β-Amyloid (35-42)****Cat. No.:** RP20145**Size:** 1 mg**Alias:** Beta-Amyloid (35-42)**Description:**

This synthetic peptide consists of amino acids 35 to 42 of b-amyloid protein. It is widely used as an immunogen for the development of the anti-b-amyloid antibodies.

**Sequence (one-letter code):**

MVGGVVIA

**Sequence (three-letter code):**

{Met}{Val}{Gly}{Gly}{Val}{Val}{Ile}{Ala}

**Solubility:** Soluble in NH<sub>4</sub>OH.**Formula:** C<sub>33</sub>H<sub>60</sub>N<sub>8</sub>O<sub>9</sub>S**Molecular Weight:** 745**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

**Note:**

\*For Non-Clinical Research Use Only \*

