

β-Amyloid (22-40)**Cat. No.:** RP20150**Size:** 1 mg**Alias:** Beta-Amyloid (22-40)**Description:**

This synthetic peptide consists of amino acids 22 to 40 of beta amyloid protein. This peptide sequence is often used in beta amyloid structure and aggregation studies.

Sequence (one-letter code):

EDVGSNKGAIIGLMVGGVV

Sequence (three-letter code):

{Glu}{Asp}{Val}{Gly}{Ser}{Asn}{Lys}{Gly}{Ala}{Ile}{Ile}{Gly}{Leu}
{Met}{Val}{Gly}{Gly}{Val}{Val}

Solubility: Soluble in NH₄OH.**Formula:** C₇₈H₁₃₅N₂₁O₂₆S**Molecular Weight:** 1,815.1**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

