

BMf-BH3**Cat. No.:** RP20201**Size:** 1 mg**Alias:** BMf-BH3**Description:**

This peptide belongs to the Bcl-2 apoptosis mediator family. Bmf is a key player in histone deacetylase (HDAC) inhibition, which alters the balance between deacetylation and acetylation.

Sequence (one-letter code):

LQHRAEVQIARKLQCIADQFHRLHT

Sequence (three-letter code):

{Leu}{Gln}{His}{Arg}{Ala}{Glu}{Val}{Gln}{Ile}{Ala}{Arg}{Lys}{Leu}
{Gln}{Cys}{Ile}{Ala}{Asp}{Gln}{Phe}{His}{Arg}{Leu}{His}{Thr}

Formula: C₁₃₁H₂₁₄N₄₅O₃₅S**Molecular Weight:** 3,012.5**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

