

DATASHEET

Version: 2016-08-18

Tetradecapeptide Renin Substrate, human**Cat. No.:** RP20216-1**Size:** 1 mg**Alias:** Preangiotensinogen (1-14) (human); Tetradecapeptide Renin Substrate (human); Angiotensinogen (1-14) (human)**Description:**

none

Sequence (one-letter code):

DRVYIHPFHLVIHN

Sequence (three-letter code): $\{\text{Asp}\}\{\text{Arg}\}\{\text{Val}\}\{\text{Tyr}\}\{\text{Ile}\}\{\text{His}\}\{\text{Pro}\}\{\text{Phe}\}\{\text{His}\}\{\text{Leu}\}\{\text{Val}\}\{\text{Ile}\}\{\text{His}\}$
 $\{\text{Asn}\}$ **Formula:** C₈₃H₁₂₂N₂₄O₁₉**Molecular Weight:** 1,760**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

