

Rev01
Update: Dec,28,2022

DATASHEET

LCMV gp33-41

Cat. No.: RP20257CN

Overview

Description	This is the H-2Db restricted epitope derived from the lymphocytic choroemeningitis virus (LCMV) glycoprotein gp 33; residues 33 to 41.
Sequence	{Lys}{Ala}{Val}{Tyr}{Asn}{Phe}{Ala}{Thr}{Met}
Sequence Shortening	KAVYNFATM
Molecular Formula	C ₄₈ H ₇₃ N ₁₁ O ₁₃ S
Molecular Weight	1044.25

Properties

Purity	> 95%
Form	Lyophilized