

DATASHEET

Version: 2016-10-10

AGA-(C8R) HNG17, Humanin derivative**Cat. No.:** RP20314**Size:** 1 mg**Description:**

This peptide is a potent humanin (HN) derivative. It completely suppresses neuronal cell death by Alzheimer's disease-relevant insults.

Sequence (one-letter code):

PAGASRLLLLTGEIDL P

Sequence (three-letter code):

{Pro}{Ala}{Gly}{Ala}{Ser}{Arg}{Leu}{Leu}{Leu}{Leu}{Thr}{Gly}{Glu}{Ile}{Asp}{Leu}{Pro}

Formula: C₇₈H₁₃₄N₂₀O₂₄**Molecular Weight:** 1,736.03**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

860 Centennial Ave., Piscataway, NJ 08854, USAToll-Free: 1-877-436-7274 Tel: 1-732-885-9188 Fax: 1-732-210-0262 Email: product@genscript.com Web: www.genscript.com

