

TIP 39, Tuberoinfundibular Neuropeptide**Cat. No.:** RP20322-1**Size:** 1 mg**Description:**

This is a tuberoinfundibular neuropeptide and parathyroid hormone 2(PTH 2)-receptor agonist from hypothalamus. Synthetic TIP39 activates human and rat PTH2 receptors.

Sequence (one-letter code):

SLALADDAAFRRERARLLAALERRHWLNSYMHKLLVLDAP

Sequence (three-letter code):

{Ser}{Leu}{Ala}{Leu}{Ala}{Asp}{Asp}{Ala}{Ala}{Phe}{Arg}{Glu}{Arg}{Ala}{Arg}{Leu}{Leu}{Ala}{Ala}{Leu}{Glu}{Arg}{Arg}{His}{Trp}{Leu}{Asn}{Ser}{Tyr}{Met}{His}{Lys}{Leu}{Leu}{Val}{Leu}{Asp}{Ala}{Pro}

Formula: C₂₀₂H₃₂₅N₆₁O₅₄S**Molecular Weight:** 4,504.2**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

