

**Protease-Activated Receptor-2, amide****Cat. No.:** RP20334-100**Size:** 100 mg**Description:**

This protease-activated receptor activating peptide (PAR-2-AP) corresponds to the PAR2 tethered ligand.

**C-Terminal:** NH<sub>2</sub>**Sequence (one-letter code):**SLIGKV-NH<sub>2</sub>**Sequence (three-letter code):**{Ser}{Leu}{Ile}{Gly}{Lys}{Val}-NH<sub>2</sub>**Formula:** C<sub>28</sub>H<sub>54</sub>N<sub>8</sub>O<sub>7</sub>**Molecular Weight:** 614.8**Storage:**

Store at -20°C.

**Note:**

\*For Non-Clinical Research Use Only \*

