

Rev01  
Update: Dec,28,2022

**DATASHEET**

# Cys-TAT(47-57)

Cat. No.: RP20343CN

## Overview

<b>Description</b>	A TAT(Transactivated-transcripton) peptide derived from the HIV TAT protein has the ability to transduce peptides or proteins into various cells.
<b>Sequence</b>	{Cys}{Tyr}{Gly}{Arg}{Lys}{Lys}{Arg}{Arg}{Gln}{Arg}{Arg}{Arg}-NH <sub>2</sub>
<b>Sequence Shortening</b>	CYGRKKRRQRRR-NH <sub>2</sub>
<b>C Terminal</b>	NH <sub>2</sub>
<b>Molecular Weight</b>	1661.99

## Properties

<b>Purity</b>	> 95%
<b>Form</b>	Lyophilized
<b>Storage</b>	Store at -20°C.