

**DATASHEET**

Version: 2016-08-18

**HEX3****Cat. No.:** RP20346**Size:** 1 mg**Description:**

This peptide is a fragment of the adenoviral hexon. Hexon is the major capsid protein of adenovirion and is comprised of three identical polypeptide chains.

**Sequence (one-letter code):**

KYSPSNVKI

**Sequence (three-letter code):**

{Lys}{Tyr}{Ser}{Pro}{Ser}{Asn}{Val}{Lys}{Ile}

**Molecular Weight:** 1,035.2**Purity:** > 95%**Storage:**

Store at -20°C.

**Note:**

\*For Non-Clinical Research Use Only \*

