

DATASHEET

Version: 2016-08-18

CDK2**Cat. No.:** RP20349**Size:** 1 mg**Description:**

This peptide is useful as a substrate for cyclin-dependent kinase 2 (cdc2) in in-vitro kinase assays.

Sequence (one-letter code):

HHASPRK

Sequence (three-letter code):

{His}{His}{Ala}{Ser}{Pro}{Arg}{Lys}

Molecular Weight: 831.93**Purity:** > 95%**Storage:**

Store at -20°C.

Note:

*For Non-Clinical Research Use Only *

