

Rev03  
Update: Feb,10,2022

## DATASHEET

# (Arg)9

Cat. No.: RP20403

## Overview

<b>Description</b>	This Arginine rich cell-penetrating peptide is capable of traversing the plasma membranes of eukaryotic cells.
<b>Sequence</b>	{Arg}{Arg}{Arg}{Arg}{Arg}{Arg}{Arg}{Arg}
<b>Sequence Shortening</b>	RRRRRRRRR
<b>Molecular Formula</b>	C <sub>54</sub> H <sub>110</sub> N <sub>36</sub> O <sub>10</sub>
<b>Molecular Weight</b>	1423.7

## Properties

<b>Purity</b>	> 95%
<b>Form</b>	Lyophilized
<b>Storage</b>	Store at -20°C. Keep tightly closed. Store in a cool dry place.