

**MAGE-A3 (195-203)****Cat. No.:** RP20413**Size:** 1 mg**Description:**

This peptide MAGE-3-encoded HLA-A24 epitope is a high MHC binder. The identification of this MAGE-3/HLA-A24 peptide may as a result effectively offer the opportunities to design peptide-based immunotherapeutic approaches that might prove to be potent in treating patients with MAGE-3-positive malignant tumors.

**Sequence (one-letter code):**

IMPKAGLLI

**Sequence (three-letter code):**

{Ile}{Met}{Pro}{Lys}{Ala}{Gly}{Leu}{Leu}{Ile}

**Formula:** C<sub>45</sub>H<sub>82</sub>N<sub>10</sub>O<sub>10</sub>S**Molecular Weight:** 955.3**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

**Note:**

\*For Non-Clinical Research Use Only \*

