

Rev03
Update: Feb,10,2022

DATASHEET

SV40 large T antigen NLS

Cat. No.: RP20414

Overview

Description	This peptide is generated from Large T antigen residue 47 to 55. It enables protein import into cell nucleus.
Sequence	{Cys}{Gly}{Gly}{Gly}{Pro}{Lys}{Lys}{Lys}{Arg}{Lys}{Val}{Glu}{Asp}
Sequence Shortening	CGGGPKKKRKVED
Molecular Formula	$C_{58}H_{104}N_{20}O_{18}S$
Molecular Weight	1401.7

Properties

Purity	> 95%
Form	Lyophilized
Storage	Store at -20°C. Keep tightly closed. Store in a cool dry place.