

DATASHEET

Version: 2016-08-18

HIV-1 TAT (48-60)**Cat. No.:** RP20423**Size:** 1 mg**Description:**

This is one of the cell-penetrating peptides generated from the human immunodeficient virus -1 Tat protein residue 48-60.

Sequence (one-letter code):

GRKKRRQRRRPPQ

Sequence (three-letter code):

{Gly}{Arg}{Lys}{Lys}{Arg}{Arg}{Gln}{Arg}{Arg}{Pro}{Pro}{Gln}

Formula: C₇₀H₁₃₁N₃₅O₁₆**Molecular Weight:** 1,719**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

860 Centennial Ave., Piscataway, NJ 08854, USA

Toll-Free: 1-877-436-7274 Tel: 1-732-885-9188 Fax: 1-732-210-0262 Email: product@genscript.com Web: www.genscript.com

