

PSA1 (141-150)**Cat. No.:** RP20434**Size:** 1 mg**Description:**

This Prostate Specific Antigen 1 peptide was used in the immunotherapy of cancer experiments.

Sequence (one-letter code):

FLTPKKLQCV

Sequence (three-letter code):

{Phe}{Leu}{Thr}{Pro}{Lys}{Lys}{Leu}{Gln}{Cys}{Val}

Formula: C₅₅H₉₃N₁₃O₁₃S**Molecular Weight:** 1,176.5**Purity:** > 95%**Storage:**

Store at -20°C. Keep tightly closed. Store in a cool dry place.

Note:

*For Non-Clinical Research Use Only *

