

Rev03
Update: Feb,10,2022

DATASHEET

[Glu1]-Fibrinopeptide B

Cat. No.: RP20439

Overview

Synonyms	Glufib
Description	This peptide is generated from fibrinopeptide B amino acid residues 1-14. It is used as a mass spec (MS) standard in proteomic research.
Sequence	{Glu}{Gly}{Val}{Asn}{Asp}{Asn}{Glu}{Glu}{Gly}{Phe}{Phe}{Ser}{Ala}{Arg}
Sequence Shortening	EGVNDNEEGFFSAR
Molecular Formula	C ₆₆ H ₉₅ N ₁₉ O ₂₆
Molecular Weight	1570.6

Properties

Purity	> 95%
Form	Lyophilized
Storage	Store the peptide at -20°C.