

Bak BH3**Cat. No.:** RP20474**Size:** 1 mg**Alias:** Bak BH3**Description:**

This peptide, generated from the BH3 domain of Bak (Flu-BakBH3), has been shown to have high-affinity binding to a surface pocket of the Bcl-XL protein that is required for its death antagonist function.

Sequence (one-letter code):

GQVGRQLAIIGDDINR

Sequence (three-letter code):

{Gly}{Gln}{Val}{Gly}{Arg}{Gln}{Leu}{Ala}{Ile}{Ile}{Gly}{Asp}{Asp}
{Ile}{Asn}{Arg}

Formula: C₇₂H₁₂₅N₂₅O₂₄**Molecular Weight:** 1,724.9**Purity:** > 95%**Storage:**

Store the peptide at -20°C.

Note:

*For Non-Clinical Research Use Only *

