

DATASHEET

Version: 2016-08-18

FliC,Serotype a (427-441) S.paratyphi A**Cat. No.:** RP20500**Size:** 1 mg**Alias:** FliC,Serotype a (427-441) S.paratyphi A**Description:**

This is amino acids 427 to 441 fragment belongs to the FliC, serotype a of the S. FliC epitope recognized by a CD4+ T cell from naturally infected C57BL/6 mice is located in a region common to a number of distinct Salmonella serovars such as the S. paratyphi. This dominant epitope is contained within the constant C-terminus of FliC.

Sequence (one-letter code):VQNRFN~~S~~AITNLGNT**Sequence (three-letter code):**

{Val}{Gln}{Asn}{Arg}{Phe}{Asn}{Ser}{Ala}{Ile}{Thr}{Asn}{Leu}{Gly}{Asn}{Thr}

Formula: C₆₉H₁₁₃N₂₃O₂₄**Molecular Weight:** 1,648.8**Purity:** > 95%**Storage:**

Store the peptide at -20°C.

Note:

*For Non-Clinical Research Use Only *

