

Rev01
Update: Dec,28,2022

DATASHEET

SV40 T-Ag-derived NLS peptide

Cat. No.: RP20513CN

Overview

Description	This peptide, a nuclear localization signal DNA tagged to this peptide efficiently translocates into the cell nucleus.
Sequence	{Pro}{Lys}{Lys}{Lys}{Arg}{Lys}{Val}{Glu}{Asp}{Pro}{Tyr}{Cys}
Sequence Shortening	PKKKRKVEDPYC
Molecular Formula	C ₆₆ H ₁₁₁ N ₁₉ O ₁₈ S ₁
Molecular Weight	1490.8

Properties

Purity	> 95%
Form	Lyophilized
Storage	Store the peptide at -20°C.