

Rev01
Update: Dec,28,2022

DATASHEET

SARS-CoV-2 ORF7b

Cat. No.: RP30015

Overview

Description	This pool includes 8 peptides derived from a peptide scan (15mers with 11 aa overlap) through the entire Non-structural Protein 7b (Protein ID: P0DTD8) of SARS-CoV-2 (Severe Acute Respiratory Syndrome-related coronavirus 2) for T cell assays.
Sequence	MIELSLIDFYLCFLAFLLLFLVLIIMLIIFWFSLELQDHNETCHA

Properties

Purity	Crude
Solubility	Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.
Form	Lyophilized
Storage	Store at -20°C.
Note	The peptides of this product are supplied as trifluoroacetate salts

Applications

T-cell assays, Immune monitoring, Antigen specific T-cell stimulation, T-cell expansion, Cellular immune response

Sequence (one-letter-code)

Code	Sequence	Code	Sequence
peptide_1	MIELSLIDFYLCFLA	peptide_2	SLIDFYLCFLAFLLF
peptide_3	FYLCFLAFLLLFLVLI	peptide_4	FLAFLLLFLVLIIMLI
peptide_5	LLFLVLIIMLIIFWFS	peptide_6	VLIIMLIIFWFSLELQ
peptide_7	LIIFWFSLELQDHNE	peptide_8	WFSLELQDHNETCHA