

Rev02
Update: Jun,15,2026

DATASHEET

SARS-CoV-2 ORF10

Cat. No.: RP30018

Overview

Description	This pool includes 7 peptides derived from a peptide scan (15mers with 11 aa overlap) through the entire ORF10 Protein (Protein ID: A0A663DJA2) of SARS-CoV-2 (Severe Acute Respiratory Syndrome-related coronavirus 2) for T cell assays (e.g. ELISPOT).
--------------------	---

Properties

Purity	Crude
Solubility	Dissolve in a minimum amount of pure DMSO (approx. 40 µl) and dilute with PBS buffer to the final concentration. Please note that the final concentration of DMSO must be below 1 % (v/v) to avoid toxicity in the biological system.
Form	Lyophilized
Storage	Store at -20°C.
Note	The peptides of this product are supplied as trifluoroacetate salts

Applications

T-cell assays, Immune monitoring, Antigen specific T-cell stimulation, T-cell expansion, Cellular immune response

Sequence (one-letter-code)

Code	Sequence	Code	Sequence
peptide_1	MGYINVFAFPFTIYS	peptide_2	NVFAFPFTIYSLLLC
peptide_3	FPFTIYSLLLCRMNS	peptide_4	IYSLLLCRMNSRNYI
peptide_5	LLCRMNSRNYIAQVD	peptide_6	MNSRNYIAQVDVWNF
peptide_7	NYIAQVDVWNFNLT		

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China