

RP30230CN
Update: Jul,13,2023

DATASHEET

SARS-CoV-2 Spike-derived Peptide (YLQPRTFLL)

Cat. No.: RP30230CN

Overview

Description	SARS-CoV-2 Spike-derived peptide S269-277 (YLQPRTFLL) is presented by HLA-A*02:01.
Sequence	{Tyr}{Leu}{Gln}{Pro}{Arg}{Thr}{Phe}{Leu}{Leu}
Sequence Shortening	YLQPRTFLL

Properties

Purity	>95%
Solubility	Can be dissolved in DMSO
Form	Lyophilized
Storage	Store at -20°C. Keep tightly closed. Store in a cool dry place.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China