

RP30231CN
Update: Jul,13,2023**DATASHEET**

SARS-CoV-2 Spike-derived Peptide (YNYKLPDDFTGCVIA)

Cat. No.: RP30231CN

Overview

Description	This peptide is derived from the entire Spike glycoprotein (Protein ID: P0DTC2) of SARS-CoV-2 (Severe Acute Respiratory Syndrome-related coronavirus 2) for T cell assays. It can be used to reveal robust pre-existing T-cell immunity in unexposed individuals.
Sequence	{Tyr}{Asn}{Tyr}{Lys}{Leu}{Pro}{Asp}{Asp}{Phe}{Thr}{Gly}{Cys}{Val}{Ile}{Ala}
Sequence Shortening	YNYKLPDDFTGCVIA

Properties

Purity	>95%
Solubility	Can be dissolved in DMSO
Form	Lyophilized
Storage	Store at -20°C. Keep tightly closed. Store in a cool dry place.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China