

Rev03
Update: Jun,09,2025

DATASHEET

SARS-CoV-2 Spike-derived Peptide (NGVEGFNCYFPLQSY)

Cat. No.: RP30234

Overview

Description	Antibody screening, T-cell assays, Immune monitoring, Antigen specific T-cell stimulation, Cellular immune response
Sequence	{Asn}{Gly}{Val}{Glu}{Gly}{Phe}{Asn}{Cys}{Tyr}{Phe}{Pro}{Leu}{Gln}{Ser}{Tyr}
Sequence Shortening	NGVEGFNCYFPLQSY
Molecular Weight	1737.89

Properties

Purity	>95%
Solubility	Soluble in dimethyl sulfoxide(Analytical grade) under 1mg/ml
Form	Lyophilized
Storage	Store at -20°C. Keep tightly closed. Store in a cool dry place.

For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China