

RP30235CN
Update: Jul,13,2023

DATASHEET

KRAS WT Peptide (VVGAGGVGK)

Cat. No.: RP30235CN

Overview

Description	KRAS peptide (VVGAGGVGK) can be used to study the crystal structure of a TCR in complex with HLA-A*11:01.
Sequence	{Val}{Val}{Val}{Gly}{Ala}{Gly}{Gly}{Val}{Gly}{Lys}
Sequence Shortening	VVGAGGVGK

Properties

Purity	>95%
Solubility	Soluble in water
Form	Lyophilized
Storage	Store at -20°C. Keep tightly closed. Store in a cool dry place.

生产商：南京金斯瑞生物科技有限公司 江苏省南京市江宁区科学园雍熙路28号

Manufacturer: Nanjing GenScript Biotech Co., Ltd. No. 28Yongxi Road, Jiangning District, Nanjing, Jiangsu, China